









## Biopsy Forceps








- Tissue sampling for alimentary tract in observation of endoscope.
- For single use only, to avoid cross infection.
- Smooth sheath to avoid any possible damage to the endoscope channel.
- The markers on PE coating can assist operator to observe the instrument place.
- Laser welding to ensure the structure strength.

# Specifications

2.3mm Biopsy Forceps	Code	Description	Tube Diameter	Length
	2316C	PE coating, w/o spike	2.3mm	1600mm
	2318C	PE coating, w/o spike	2.3mm	1800mm
	2323C	PE coating, w/o spike	2.3mm	2300mm
	2316	w/o spike	2.3mm	1600mm
	2318	w/o spike	2.3mm	1800mm
	2323	w/o spike	2.3mm	2300mm
	2316SC	PE coating, with spike	2.3mm	1600mm
	2318SC	PE coating, with spike	2.3mm	1800mm
	2323SC	PE coating, with spike	2.3mm	2300mm
	2316S	with spike	2.3mm	1600mm
	2318S	with spike	2.3mm	1800mm
	2323S	with spike	2.3mm	2300mm
	2316TC	PE coating, alligator, w/o spike	2.3mm	1600mm
	2318TC	PE coating, alligator, w/o spike	2.3mm	1800mm
	2323TC	PE coating, alligator, w/o spike	2.3mm	2300mm
	2316T	alligator, w/o spike	2.3mm	1600mm
	2318T	alligator, w/o spike	2.3mm	1800mm
	2323T	alligator, w/o spike	2.3mm	2300mm
	2316TSC	PE coating, alligator, with spike	2.3mm	1600mm
	2318TSC	PE coating, alligator, with spike	2.3mm	1800mm
	2323TSC	PE coating, alligator, with spike	2.3mm	2300mm
	2316TS	alligator, with spike	2.3mm	1600mm
	2318TS	alligator, with spike	2.3mm	1800mm
	2323TS	alligator, with spike	2.3mm	2300mm

# Specifications

1.8mm Biopsy Forceps	Code	Description	Tube Diameter	Length
	1812C	PE coating, w/o spike	1.8mm	1200mm
	1816C	PE coating, w/o spike	1.8mm	1600mm
	1818C	PE coating, w/o spike	1.8mm	1800mm
	1812	w/o spike	1.8mm	1200mm
	1816	w/o spike	1.8mm	1600mm
	1818	w/o spike	1.8mm	1800mm
	1812TC	PE coating, alligator, w/o spike	1.8mm	1200mm
	1816TC	PE coating, alligator, w/o spike	1.8mm	1600mm
	1818TC	PE coating, alligator, w/o spike	1.8mm	1800mm
	1812T	alligator, w/o spike	1.8mm	1200mm
	1816T	alligator, w/o spike	1.8mm	1600mm
	1818T	alligator, w/o spike	1.8mm	1800mm
	1812SC	PE coating, with spike	1.8mm	1200mm
	1816SC	PE coating, with spike	1.8mm	1600mm
	1818SC	PE coating, with spike	1.8mm	1800mm

Jumbo Biopsy Forceps{3.0mm}	Code	Description	Tube Diameter	Length
	3018C	PE coating, w/o spike	3.0mm	1800mm
	3023C	PE coating, w/o spike	3.0mm	2300mm
	3018	w/o spike	3.0mm	1800mm
	3023	w/o spike	3.0mm	2300mm

\* Marker type is available for all PE coating biopsy forceps.

# Cleaning Brushes



Specifications	Code	Description	Tube Diameter	Length	Use
	CBV10051723B	5.0mm bristle+10mm bristle with plastic ball with plastic handle	1.7mm	2300mm	Endoscope Valve Gastroscope/Colonoscope
	CB100000	10mm bristle plastic handle	-	-	Endoscope Valve
	CB051000	5.0mm bristle+10mm bristle, plastic handle	-	-	Endoscope Valve
	CB02001010B	2.0mm bristle x 1 with plastic ball	1.0mm	1000mm	Bronchoscope
	CB03001718B	3.0mm bristle x 1 with plastic ball	1.7mm	1800mm	Bronchoscope
	CB05001723B	5.0mm bristle x 1 with plastic ball	1.7mm	2300mm	Gastroscope/Colonoscope
	CB02021018B	2.0mm bristle x 2 with plastic ball	1.0mm	1000mm	Bronchoscope
	CB03031718B	3.0mm bristle x 2 with plastic ball	1.7mm	1800mm	Bronchoscope
	CB05051723B	5.0mm bristle x 2 with plastic ball	1.7mm	2300mm	Gastroscope/Colonoscope
	CB05101723B	5.0mm bristle+10mm bristle with plastic ball	1.7mm	2300mm	Endoscope Valve Gastroscope/Colonoscope